



New and In Stock Now!
The next level of
True Blue Performance

Topaz GC inlet liners feature revolutionary technology and inertness to deliver the next level of:

- Deactivation
- Reproducibility
- Productivity

www.restek.com/topaz



* For details, see back page.



RESTEK

Pure Chromatography















www.restek.com/topaz

The next level of True Blue Performance

Topaz GC inlet liners feature revolutionary technology and inertness to deliver the next level of True Blue Performance:

- **Deactivation**—unbelievably low breakdown for accurate and precise low-level GC analyses.
- **Reproducibility**—unbeatable manufacturing controls and QC testing for superior reliability across compound classes.
- **Productivity**—unparalleled cleanliness for maximized GC uptime and lab throughput.
- **100% Satisfaction**—if a liner doesn't perform to your expectations, we will replace it or credit your account.

Topaz Inlet Liners for Agilent GCs

Split Liners for Agilent GCs	ID OD x Length	Similar to Agilent Part #	qty.	cat.#
 RESTEK 1 mm Split	1.0 mm 6.3 mm x 78.5 mm		5-pk.	23333
 RESTEK 4 mm Straight w/Wool	4.0 mm 6.3 mm x 78.5 mm	5183-4691, 5190-3164	5-pk.	23300
 RESTEK 2 mm Precision Liner w/Wool	2.0 mm 6.3 mm x 78.5 mm		5-pk.	23468
 RESTEK 4 mm Precision Liner w/Wool	4.0 mm 6.3 mm x 78.5 mm	210-4004-5	5-pk.	23305
 RESTEK 4 mm Cyclo	4.0 mm 6.3 mm x 78.5 mm		5-pk.	23312
Splitless Liners for Agilent GCs	ID OD x Length	Similar to Agilent Part #	qty.	cat.#
 RESTEK 2 mm Splitless	2.0 mm 6.5 mm x 78.5 mm	5183-4703	5-pk.	23313
 RESTEK 2 mm Splitless w/Wool	2.0 mm 6.5 mm x 78.5 mm		5-pk.	23314
 RESTEK 2 mm Single Taper	2.0 mm 6.5 mm x 78.5 mm		5-pk.	23315
 RESTEK 2 mm Single Taper w/Wool	2.0 mm 6.5 mm x 78.5 mm		5-pk.	23316
 RESTEK 4 mm Straight	4.0 mm 6.3 mm x 78.5 mm	210-3003-05	5-pk.	23301
 RESTEK 4 mm Straight w/Wool	4.0 mm 6.3 mm x 78.5 mm	5183-4691, 5190-3164	5-pk.	23300
 RESTEK 4 mm Single Taper	4.0 mm 6.5 mm x 78.5 mm	5183-4695	5-pk.	23302
 RESTEK 4 mm Single Taper w/Wool	4.0 mm 6.5 mm x 78.5 mm	5183-4693, 5190-3163	5-pk.	23303
 RESTEK 4 mm Double Taper	4.0 mm 6.5 mm x 78.5 mm	5183-4705	5-pk.	23308
 RESTEK 4 mm Cyclo Double Taper	4.0 mm 6.5 mm x 78.5 mm		5-pk.	23310

Patented















COLUMN INSTALLS THIS END

Topaz inlet liners feature a 100% satisfaction guarantee.

All Topaz inlet liners are produced with proprietary Restek processes and apply some of the tightest testing specifications in the industry. The exceptionally inert Topaz liners analysts want most are in stock now and are guaranteed to deliver outstanding performance that exceeds your expectations.



Topaz Inlet Liners for Agilent GCs

Split/Splitless Liners for Agilent GCs	ID OD x Length	Similar to Agilent Part #	qty.	cat.#
 4 mm Straight w/Wool	4.0 mm 6.3 mm x 78.5 mm	5183-4691, 5190-3164	5-pk.	23300
 Low Pressure Drop Precision Liner w/Wool	4.0 mm 6.3 mm x 78.5 mm		5-pk.	23309
 Low Pressure Drop Liner w/Wool	4.0 mm 6.3 mm x 78.5 mm	5183-4701, 5183-4712, 5190-3165	5-pk.	23467
CIS4 and PTV Liners for Agilent GCs	ID OD x Length	Similar to Agilent Part #	qty.	cat.#
 On Column PTV	1.7 mm 3.0 mm x 71 mm		5-pk.	23430
 Single Baffle PTV	2.0 mm 3.0 mm x 71 mm	5183-2036	10-pk.	23431
 Single Baffle PTV w/Wool	2.0 mm 3.0 mm x 71 mm	5183-2038	10-pk.	23432
 Baffled PTV	1.5 mm 3.0 mm x 71 mm	5183-2037	10-pk.	23433
CIS4/TDU Liners for Agilent GCs	ID OD x Length	Similar to Gerstel Part #	qty.	cat.#
 Single Baffle CIS4/TDU	2.0 mm 3.0 mm x 78 mm	013775-010-00	10-pk.	23283
 Single Baffle CIS4/TDU w/Wool	2.0 mm 3.0 mm x 78 mm	012742-010-00	10-pk.	23284
 Baffled CIS4/TDU	2.0 mm 3.0 mm x 78 mm	012436-010-00	10-pk.	23282
Analytical Controls TPI Inlet Liners for Agilent GCs	ID OD x Length		qty.	cat.#
 TPI for 0.53 mm ID columns	2.4 mm 4.0 mm x 71 mm		10-pk.	23429
 TPI for 0.25/0.32 mm ID columns	2.4 mm 4.0 mm x 71 mm		10-pk.	23428
SPME Liners for Agilent GCs	ID OD x Length		qty.	cat.#
 SPME Liner	0.75 mm 6.35 mm x 78.5 mm		5-pk.	23434
 SPME Liner	1.8 mm 6.5 mm x 78.5 mm		5-pk.	23280

Patented

COLUMN INSTALLS THIS END

QUANTITY DISCOUNTS AVAILABLE for orders of five or more 5-packs!

Topaz Inlet Liners for Agilent GCs

Direct Injection Liners for Agilent GCs (for 0.25/0.32/0.53 mm ID Columns)	ID OD x Length	qty.	cat.#
RESTEK Drilled Uniliner (hole near bottom)	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23306
RESTEK Drilled Uniliner (hole near bottom) w/Wool	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23307
RESTEK Drilled Uniliner (hole near top)	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23311

Topaz Inlet Liners for Bruker/Varian GCs

Liners for Bruker/Varian 1177 S/SL Injection Ports	ID OD x Length	qty.	cat.#
RESTEK 4 mm Split Liner w/Glass Frit	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23330
RESTEK 4 mm Precision Liner w/Wool	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23328
RESTEK 4 mm Single Taper	4.0 mm 6.5 mm x 78.5 mm	5-pk.	23331
RESTEK 4 mm Single Taper w/Wool	4.0 mm 6.5 mm x 78.5 mm	5-pk.	23332

SPI Liners for Bruker/Varian S/SL Injection Ports	ID OD x Length	Similar to Bruker/Varian Part #	qty.	cat.#
RESTEK SPI Liner	0.53 mm 4.6 mm x 54 mm		5-pk.	23460
RESTEK SPI Liner	0.80 mm 4.6 mm x 54 mm	4984, 5122991	5-pk.	23461

Liners for Bruker/Varian 1078/1079 Injection Ports	ID OD x Length	qty.	cat.#
RESTEK 3.4 mm Split-No Frit	3.4 mm 5.0 mm x 54 mm	5-pk.	23329
RESTEK Split w/Glass Frit	3.4 mm 5.0 mm x 54 mm	5-pk.	23462
RESTEK Splitless	2.0 mm 5.0 mm x 54 mm	5-pk.	23463
RESTEK SPME Liner	0.75 mm 5.0 mm x 54 mm	5-pk.	23465
RESTEK Split Precision Liner w/Wool	3.4 mm 5.0 mm x 54 mm	5-pk.	23466

Patented

COLUMN INSTALLS THIS END

















In stock now — order today at www.restek.com/topaz


Topaz Inlet Liners for DANI GCs



Liners for DANI GCs	ID OD x Length	qty.	cat.#
 RESTEK	0.8 mm		
0.8 mm ID Straight Liner	6.0 mm x 70 mm	5-pk.	23293
0.8 mm ID Straight Liner (includes: ferrule set)		5-pk.	23286
 RESTEK	2.0 mm		
2.0 mm ID Straight Liner	6.0 mm x 70 mm	5-pk.	23292
2.0 mm ID Straight Liner (includes: ferrule set)		5-pk.	23285
 RESTEK	2.2 mm		
2.2 mm ID Straight Liner w/Wool	4.0 mm x 85.7 mm	5-pk.	23291
 RESTEK	2.2 mm		
2.2 mm ID Dimpled Liner w/Wool	4.0 mm x 55 mm	5-pk.	23296
2.2 mm ID Dimpled Liner w/Wool (includes: ferrule set)		5-pk.	23289
 RESTEK	4.0 mm		
4.0 mm ID Straight Liner	6.0 mm x 70 mm	5-pk.	23294
4.0 mm ID Straight Liner (includes: ferrule set)		5-pk.	23287
 RESTEK	4.0 mm		
4.0 mm ID Straight Liner w/Wool	6.0 mm x 70 mm	5-pk.	23295
4.0 mm ID Straight Liner w/Wool (includes: ferrule set)		5-pk.	23288
 RESTEK	4.0 mm		
4.0 mm ID Single Taper Liner w/Wool	6.0 mm x 70 mm	5-pk.	23297
4.0 mm ID Single Taper Liner w/Wool (includes: ferrule set)		5-pk.	23290

Topaz Inlet Liners for PerkinElmer GCs

Split Liners for PerkinElmer GCs	ID OD x Length	qty.	cat.#
 RESTEK	4.0 mm		
Splitter w/Wool	6.2 mm x 92.1 mm	5-pk.	23449
 RESTEK	4.0 mm		
4 mm Split Precision Liner w/Wool	6.2 mm x 92.1 mm	5-pk.	23450

Splitless Liners for PerkinElmer GCs	ID OD x Length	qty.	cat.#
 RESTEK	2.0 mm		
2 mm Splitless w/Wool	6.2 mm x 92.1 mm	5-pk.	23451

PSS Liners for PerkinElmer GCs	ID OD x Length	qty.	cat.#
 RESTEK	2.0 mm		
Auto SYS XL PSS Split/Splitless w/Wool	4.0 mm x 86.2 mm	5-pk.	23317






DI Liners for PerkinElmer GCs	ID OD x Length	qty.	cat.#
 RESTEK	4.0 mm		
Open-Top Uniliner w/Wool	6.2 mm x 92.1 mm	5-pk.	23452
 RESTEK	2.0 mm		
PSS Drilled Uniliner (hole near top)	4.0 mm x 86.2 mm	5-pk.	23281
 RESTEK	4.0 mm		
Drilled Uniliner (hole near top)	6.2 mm x 92.1 mm	5-pk.	23453
 RESTEK	4.0 mm		
Drilled Uniliner (hole near bottom)	6.2 mm x 92.1 mm	5-pk.	23454





Patented


QUANTITY DISCOUNTS AVAILABLE for orders of five or more 5-packs!


Topaz Inlet Liners for Shimadzu GCs

COLUMN INSTALLS THIS END


Split Liners for Shimadzu 17A, 2010, and 2014 GCs	ID OD x Length	qty.	cat.#
 3.5 mm Split	3.5 mm 5.0 mm x 95 mm	5-pk.	23318
 3.5 mm Split w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23319
 3.5 mm Precision Liner w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23320
 Single Taper Precision Liner w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23455
 Cyclo	3.5 mm 5.0 mm x 95 mm	5-pk.	23456



Splitless Liners for Shimadzu 17A, 2010, and 2014 GCs	ID OD x Length	qty.	cat.#
 3.5 mm Single Taper	3.5 mm 5.0 mm x 95 mm	5-pk.	23321
 3.5 mm Single Taper w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23336
 3.5 mm Single Taper w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23322
 Double Taper	3.5 mm 5.0 mm x 95 mm	5-pk.	23457

Split/Splitless Liners for Shimadzu 17A, 2010, and 2014 GCs	ID OD x Length	qty.	cat.#
 Split/Splitless w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23458

Liners for Shimadzu 17A PTV GCs	ID OD x Length	qty.	cat.#
 PTV w/Wool	1.6 mm 4.0 mm x 95 mm	5-pk.	23435

SPME Liners for Shimadzu 17A, 2010, and 2014 GCs	ID OD x Length	qty.	cat.#
 SPME Liner	1.8 mm 5.0 mm x 95 mm	5-pk.	23279






DI Liners for Shimadzu 17A, 2010, and 2014 GCs	ID OD x Length	qty.	cat.#
 Uniliner w/Wool	3.5 mm 5.0 mm x 95 mm	5-pk.	23459






Liners for Shimadzu 2010 PTV GCs	ID OD x Length	qty.	cat.#
 PTV 2010	1.5 mm 3.5 mm x 95 mm	10-pk.	23471
 PTV 2010 w/Wool	1.5 mm 3.5 mm x 95 mm	10-pk.	23472


Patented

Topaz Inlet Liners for Thermo Scientific GCs

COLUMN INSTALLS THIS END

Split Liners for Thermo TRACE, 8000, 8000 TOP, and Focus SSL	ID OD x Length	qty.	cat.#
 Split Straight	3.0 mm 8.0 mm x 105 mm	5-pk.	23439
 Split Straight w/Wool	3.0 mm 8.0 mm x 105 mm	5-pk.	23440
 5 mm Straight	5.0 mm 8.0 mm x 105 mm	5-pk.	23323
 5 mm Straight w/Wool	5.0 mm 8.0 mm x 105 mm	5-pk.	23324
 5 mm Precision Liner w/Wool	5.0 mm 8.0 mm x 105 mm	5-pk.	23327

Splitless Liners for Thermo TRACE, 8000, 8000 TOP, and Focus SSL	ID OD x Length	qty.	cat.#
 Splitless	3.0 mm 8.0 mm x 105 mm	5-pk.	23441
 Splitless w/Wool	3.0 mm 8.0 mm x 105 mm	5-pk.	23442
 Splitless Precision Liner w/Wool	5.0 mm 8.0 mm x 105 mm	5-pk.	23443
 5 mm Splitless	5.0 mm 8.0 mm x 105 mm	5-pk.	23325
 5 mm Splitless w/Wool	5.0 mm 8.0 mm x 105 mm	5-pk.	23326

Split Liners for Thermo TRACE 1300/1310	ID OD x Length	qty.	cat.#
 Split	1.0 mm 6.3 mm x 78.5 mm	5-pk.	23448

Patented

QUANTITY DISCOUNTS AVAILABLE for orders of five or more 5-packs!



Restek Enviro

YOUR GUIDE FOR Environmental Sampling and Analysis









*LC and GC Supplies, Technical
Expertise, and Plus 1 Service*

www.restek.com/enviro



COLUMN INSTALLS THIS END

Topaz Inlet Liners for Thermo Scientific GCs

Splitless Liners for Thermo TRACE 1300/1310		ID	qty.	cat.#
		OD x Length		
	Splitless Straight	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23445
	4 mm Splitless Single Taper	4.0 mm 6.5 mm x 78.5 mm	5-pk.	23446
	4 mm Splitless Single Taper w/Wool	4.0 mm 6.5 mm x 78.5 mm	5-pk.	23447
Split/Splitless Liners for Thermo TRACE 1300/1310		ID	qty.	cat.#
		OD x Length		
	4 mm Split/Splitless Straight w/Wool	4.0 mm 6.3 mm x 78.5 mm	5-pk.	23444
Split Liners for Thermo Scientific TRACE PTV		ID	qty.	cat.#
		OD x Length		
	Split PTV	1.0 mm 2.75 mm x 120 mm	5-pk.	23436
	Split PTV	2.0 mm 2.75 mm x 120 mm	5-pk.	23437
	Baffled PTV	2.0 mm 2.75 mm x 120 mm	5-pk.	23438
SPME Liner for Thermo TRACE 1300/1310		ID	qty.	cat.#
		OD x Length		
	SPME Liner	1.8 mm 6.5 mm x 78.5 mm	5-pk.	23278

Patented



Perfect Installation Every Time

There is never a question about which way to install a Topaz liner—simply orient it so that the column installs toward the "R" in the Restek logo.

*100% SATISFACTION GUARANTEE: If your Topaz inlet liner does not perform to your expectations for any reason, simply contact Restek Technical Service or your local Restek representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty

RESTEK
Pure Chromatography

Questions about this or any other Restek product?

Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2017 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# GNSS2635-UNV